

Table 1 Results of EMR+ ex vivo * time from injection to complete lifting of lesion and flushing with warm water ** time from snare opening to extraction of resection specimen (all resections were performed in the area of the greater curvature)

Resection	Specimen L x B x H (mm)	Time injection* (min)	Time resection** (min)	Time EMR+ (min)
1	30 x 27 x 10	4	7	11
2	25 x 24 x 10.4	3	3	6
3	30 x 24 x 9.2	3	3	6
4	22 x 20 x 7	4	3	7
5	31 x 25 x 13.3	3	3	6
6	35 x 25 x 11.6	2	4	6
7	34 x 29 x 12.7	2	3	5
8	40 x 33 x 14.4	3	2	5
9	29 x 25 x 8.8	3	4	7
10	32 x 26 x 10.3	4	3	7
11	30 x 25 x 11.4	3	2	5
12	30 x 26 x 9.6	2	9	11
13	30 x 28 x 12,4	2	4	6
14	25 x 22 x 8.5	3	4	7
15	30 x 29 x 10.5	5	3	8
16	36 x 30 x 13.5	3	4	7
17	33 x 30 x 16.6	4	5	9
18	32 x 28 x 14.4	2	5	7
19	25 x 25 x 12.7	3	4	7
20	26 x 22 x 6.5	4	5	9
21	29 x 29 x 11.5	4	4	8
22	37 x 30 x 10.1	2	4	6
Median	30 x 26 x 11	3	4	7

Table 2 Results of EMR+ in vivo * time from injection to complete lifting of lesion ** time from snare opening to extraction of resection specimen (all resections were performed in the area of the greater curvature) *** management by coagulation

Resection	Specimen L x B x H (mm)	Time injection* (min)	Time resection** (min)	Time EMR+ (min)	Adverse events	Macroscopic complete
1	40 x 35 x 10	2	6	8	-	yes
2	40 x 35 x 11	6	3	9	-	yes
3	40 x 38 x 13	2	6	8	-	yes
4	40 x 38 x 13	2	10	12	-	yes
5	30 x 30 x 12	2	6	8	-	no
6	35 x 35 x 11	2	4	6	-	yes
7	40 x 40 x 11	2	2	4	-	yes
8	35 x 35 x 9	1	2	3	-	yes
9	35 x 35 x 13	2	2	4	-	yes
10	35 x 38 x 10	3	2	5	-	yes
11	40 x 35 x 14	1	2	3	-	yes
12	35 x 35 x 9	2	3	5	minor bleeding***	yes
13	35 x 35 x 8	2	2	4	-	yes
Median	35 x 35 x 11	2	3	5	-	-