WÙ

# World Journal of Gastroenterology

Submit a Manuscript: https://www.f6publishing.com

World J Gastroenterol 2024 March 7; 30(9): 1121-1131

ISSN 1007-9327 (print) ISSN 2219-2840 (online)

DOI: 10.3748/wjg.v30.i9.1121

ORIGINAL ARTICLE

## **Clinical Trials Study** Detachable string magnetically controlled capsule endoscopy for the noninvasive diagnosis of esophageal diseases: A prospective, blind clinical study

Yan-Ling Yang, Huang-Wen Qin, Zhao-Yu Chen, Hui-Ning Fan, Yi Yu, Wei Da, Jin-Shui Zhu, Jing Zhang

Specialty type: Gastroenterology and hepatology

Provenance and peer review: Unsolicited article; Externally peer

reviewed

Peer-review model: Single blind

#### Peer-review report's scientific quality classification

Grade A (Excellent): 0 Grade B (Very good): B Grade C (Good): 0 Grade D (Fair): 0 Grade E (Poor): 0

P-Reviewer: Oh SH, South Korea

Received: November 22, 2023 Peer-review started: November 22, 2023 First decision: December 7, 2023 Revised: January 14, 2024 Accepted: February 1, 2024 Article in press: February 1, 2024 Published online: March 7, 2024



Yan-Ling Yang, Huang-Wen Qin, Zhao-Yu Chen, Hui-Ning Fan, Yi Yu, Wei Da, Jin-Shui Zhu, Jing Zhang, Department of Gastroenterology, Shanghai Sixth People's Hospital Affiliated to Shanghai Jiao Tong University School of Medicine, Shanghai 200233, China

Corresponding author: Jing Zhang, MD, Associate Chief Physician, Department of Gastroenterology, Shanghai Sixth People's Hospital Affiliated to Shanghai Jiao Tong University School of Medicine, No. 600 Yishan Road, Shanghai 200233, China. jing5522724@vip.163.com

## Abstract

## BACKGROUND

Traditional esophagogastroduodenoscopy (EGD), an invasive examination method, can cause discomfort and pain in patients. In contrast, magnetically controlled capsule endoscopy (MCE), a noninvasive method, is being applied for the detection of stomach and small intestinal diseases, but its application in treating esophageal diseases is not widespread.

## AIM

To evaluate the safety and efficacy of detachable string MCE (ds-MCE) for the diagnosis of esophageal diseases.

## **METHODS**

Fifty patients who had been diagnosed with esophageal diseases were prospectively recruited for this clinical study and underwent ds-MCE and conventional EGD. The primary endpoints included the sensitivity, specificity, positive predictive value, negative predictive value, and diagnostic accuracy of ds-MCE for patients with esophageal diseases. The secondary endpoints consisted of visualizing the esophageal and dentate lines, as well as the subjects' tolerance of the procedure.

## RESULTS

Using EGD as the gold standard, the sensitivity, specificity, positive predictive value, negative predictive value, and diagnostic accuracy of ds-MCE for esophageal disease detection were 85.71%, 86.21%, 81.82%, 89.29%, and 86%, respectively. ds-MCE was more comfortable and convenient than EGD was, with 80% of patients feeling that ds-MCE examination was very comfortable or



WJG https://www.wjgnet.com

comfortable and 50% of patients believing that detachable string v examination was very convenient.

#### **CONCLUSION**

This study revealed that ds-MCE has the same diagnostic effects as traditional EGD for esophageal diseases and is more comfortable and convenient than EGD, providing a novel noninvasive method for treating esophageal diseases.

**Key Words:** Clinical trial; Detachable string magnetically controlled capsule endoscopy; Esophagogastroduodenoscopy; Noninvasive diagnosis; Esophageal diseases

©The Author(s) 2024. Published by Baishideng Publishing Group Inc. All rights reserved.

Core Tip: Our study highlights detachable string magnetically controlled capsule endoscopy (ds-MCE) as an innovative, noninvasive technique for esophageal diagnosis that matches the 86% accuracy of traditional esophagogastroduodenoscopy while significantly enhancing patient comfort. With 80% of patients reporting a favorable experience, ds-MCE stands to increase compliance and revolutionize the approach to esophageal disease management, offering a promising leap forward in gastroenterological diagnostics with minimal patient discomfort.

Citation: Yang YL, Qin HW, Chen ZY, Fan HN, Yu Y, Da W, Zhu JS, Zhang J. Detachable string magnetically controlled capsule endoscopy for the noninvasive diagnosis of esophageal diseases: A prospective, blind clinical study. World J Gastroenterol 2024; 30(9): 1121-1131

URL: https://www.wjgnet.com/1007-9327/full/v30/i9/1121.htm DOI: https://dx.doi.org/10.3748/wjg.v30.i9.1121

## INTRODUCTION

Esophageal diseases, including gastroesophageal reflux disease, Barrett's esophagus, and esophageal tumors, pose serious threats to human health[1]. Long-term chronic inflammation is strongly associated with an increased risk of Barrett's esophagus and esophageal cancer[2]. Most malignant cases can metastasize to distant sites, leading to poor prognosis[3,4]. Consequently, early diagnosis of esophageal disease is crucial for preventing the occurrence of esophageal cancer.

Esophagogastroduodenoscopy (EGD) is the "gold standard" for detecting upper gastrointestinal diseases; this method allows visualization of the mucosa and sampling, but EGD requires intubation and even anesthesia, resulting in poor compliance. Esophageal capsule endoscopy (ECE) is a well-tolerated and safe procedure[5] that was initially used to diagnose small intestinal diseases[6]. In 2004, the first-generation ECE was authorized, and video images were acquired from both ends of the device at a rate of 2 frames/second/end[7]. Second-generation ECE was approved in 2007 to provide improved image quality at 15 frames/second and high spatial resolution[8]. A third-generation ECE features a 174° angle of view and a 35 frames/second rate[9] but fails to improve detection rates[10]. ECE, as a screening tool for esophageal diseases, has shown a moderate sensitivity and specificity of 77% and 86%, respectively, of cases, whereas EGD has a sensitivity and specificity of 76% and 90% [11], respectively; however, ECE has a high false-positive rate[8], making it ineffective as an alternative to EGD. In 2013, a magnetically controlled capsule endoscopy (MCE) system was developed and approved in China<sup>[12]</sup>; this system involves a capsule manipulated by an external magnetic field to visualize the mucosa. MCE is nearly equivalent to conventional EGD due to its diagnostic accuracy<sup>[13]</sup> and high sensitivity for detecting minor erosions<sup>[12]</sup>. Another study demonstrated that the diagnostic accuracy of MCE for esophageal lesions can reach 86% [14]. Moreover, detachable string MCE (ds-MCE) was developed to control the movement of the capsule, allowing direct and repeated observation of the esophagus[15].

Our previous study showed that MCE is a safe and noninvasive endoscopic imaging method with a highly accurate detection rate for gastric and small intestinal diseases [16]. In the present study, ds-MCE was applied for the detection of esophageal diseases, and we found that ds-MCE had the same diagnostic effects as traditional EGD for esophageal diseases and was more comfortable and convenient than EGD was, providing a novel noninvasive method for detecting esophageal diseases.

#### MATERIALS AND METHODS

#### Study design

This was a prospective, blinded, self-controlled clinical study conducted at Shanghai Sixth People's Hospital from November 29th, 2019, to September 20th, 2021. This study was approved by the Ethics Committee of Shanghai Sixth People's Hospital, No. 2019-082-(1).



#### Study participants

Participants aged 18-75 years with previously diagnosed esophageal diseases were eligible for this clinical study. Patients with one of the following conditions were excluded: (1) No surgical condition or refusal to undergo any abdominal surgery; (2) an intracorporeal pacemaker; (3) intracorporeal implantation of electronic devices such as cochlear implants, magnetic metal drug infusion pumps, neurostimulators, or magnetic metal foreign bodies; (4) pregnancy; (5) dysphagia, known or suspected gastrointestinal obstruction, stricture, or fistula; significant gastrointestinal bleeding; a history of gastrointestinal surgery or history of abdominal surgery with altered gastrointestinal anatomy; or history of abdominal radiation; (6) contraindication to electrogastroscopy; and (7) other conditions that the investigator believed were not appropriate for this study. All the subjects understood and agreed to the study protocol and voluntarily signed the informed consent form.

#### Procedures

After providing informed consent, each subject underwent ds-MCE followed by EGD, which is the standard diagnostic method for comparison with ds-MCE. The EGD operating and diagnosing physicians were unblinded to the ds-MCE operating and reviewing physicians. The procedure for this study is shown in Figure 1.

The ds-MCE (specification: NU-1/CEE-1), which consisted of an MCE system and a single-use detachable string (Figure 2), was obtained from Shanghai Ankon Medical Technology Co., Ltd. (Shanghai, China). The examiner used only one capsule to complete the upper gastrointestinal examination. The detachable string of the single-use capsule endoscope was connected to one end of the capsule, and the other end was attached to a syringe, which moved the capsule up and down to examine the esophagus. Following the examination of the esophagus, air was injected through the syringe to separate the capsule from the string, and the capsule was then inserted into the stomach to examine the stomach and small intestine.

#### Study outcomes and definition

Esophageal lesions were detected by ds-MCE and EGD, and EGD was used as the gold standard for the diagnosis of esophageal disease. The primary endpoints of this study included the sensitivity, specificity, positive predictive value, negative predictive value, and diagnostic accuracy of ds-MCE for detecting esophageal disease (Table 1): (1) Sensitivity: The probability of a positive EGD and a positive MCE result, which was an indicator of the ability to detect a positive result in this test; (2) Specificity: The specificity is the probability of a negative EGD and a negative MCE, which is an indicator of the ability to detect a negative result; (3) Predictive value: Includes positive and negative predictive values. The higher the value was, the more likely the probability of the subject having the disease was to be estimated; and (4) Diagnostic accuracy: The probability that the ds-MCE was consistent with the EGD findings.

Secondary endpoints: (1) Esophageal visualization: The ability to clearly capture and preserve images of the upper, middle, and lower esophagus; (2) Dentate line visualization: Clear capture of the dentate line images; and (3) Subject tolerance: Participants completed a questionnaire to assess their level of comfort with EGD and ds-MCE, with ratings ranging from 0 to 4 (0 = very comfortable; 1 = comfortable; 2 = tolerable; 3 = uncomfortable; 4 = very uncomfortable); the ease of inspection ranged from 0 to 3 (0 = very convenient; 1 = convenient; 2 = inconvenient; 3 = very inconvenient). This study evaluated the safety of the interventions by evaluating the incidence of adverse events and complications during and after capsule endoscopy and EGD.

#### Statistical analysis

The statistical analysis of the data was conducted with SPSS 26.0 and GraphPad Prism 9. All hypothesis tests were performed at a two-sided 5% level of significance. Continuous variables were summarized using descriptive statistics. For normally distributed variables, the mean and standard deviation were calculated; for variables with a skewed distribution, the median and interquartile spacing were calculated. The Kruskal-Wallis test was used for comparisons between groups.

## RESULTS

#### Patient characteristics

A total of 50 patients (average age of 53.26 years; 46% female) were included; 48 patients had a history of esophagitis, 1 patient had a history of Barrett's esophagus, 2 patients had a history of esophageal hiatal hernia, 8 patients had superficial gastritis, and 2 patients had atrophic gastritis. The demographic data of the 50 subjects are presented in Table 2.

#### Diagnostic accuracy and safety of ds-MCE

Fifty participants successfully underwent ds-MCE and did not experience any adverse reactions or device defects associated with ds-MCE or EGD. A total of 21 patients (42%) with esophageal diseases were detected by EGD, and 22 patients (44%) were detected by ds-MCE. Among them, 3 patients were diagnosed with esophageal disease by EGD but were negative by ds-MCE, while 4 patients were diagnosed with esophageal disease by ds-MCE but were negative by EGD. The detection of specific esophageal diseases is shown in Supplementary Table 1.

Using EGD as the gold standard, the sensitivity of ds-MCE for detecting esophageal disease was 85.71% (95%CI: 62.64%, 96.24%), the specificity was 86.21% (95%CI: 67.43%, 95.49%), the positive predictive value was 81.82% (95%CI: 58.99%, 94.01%), the negative predictive value was 89.29% (95%CI: 70.63%, 97.19%), and the diagnostic concordance was

#### Yang YL et al. A prospective, blind clinical study

Table 1 Methods for detecting diseases in experimental and control groups					
		EGD			
		Positive	Negative	Total	
ds-MCE	Positive	a	b	a + b (R1)	
	Negative	c	d	c + d (R2)	
	Total	a + c (C1)	b + d (C2)	a + b + c + d(N)	

(1) Sensitivity = Number of true positives/(number of true positives + number of false negatives) × 100 = a/C1; (2) Specificity = number of true negatives/(number of true negatives + number of false positives) × 100 = d/C2; (3) Positive predictive value = number of true positives/(number of true positives + number of false positives) × 100 = a/R1; (4) Negative predictive value = number of true negatives/(number of true negatives + number of false negatives) × 100 = d/R2; And (5) Diagnostic accuracy (concordance) = (number of true positives + number of true negatives)/total number × 100 = (a + d)/ N. Based on the definitions, two-sided 95%CI were calculated for each of the above indicators. EGD: Esophagogastroduodenoscopy; ds-MCE: Detachable string magnetically controlled capsule endoscopy.




Figure 1 The procedure for this study. ds-MCE: Detachable string magnetically controlled capsule endoscopy.




Figure 2 The operational procedure of detachable string magnetically controlled capsule endoscopy. A: A single-use capsule endoscope and detachable string; B: The detachable string is connected to one end of the capsule and the other end is attached to a syringe; C: Air is injected into the hollow string through the syringe; D: Separating the capsule from the string.

Baishidena® WJG | https://www.wjgnet.com

Table 2 Demographic and baseline characteristics	
	Subjects ( <i>n</i> = 50)
Sex, % (n)	
Male	54% (27/50)
Female	46% (23/50)
Age (yr)	
n	50
Mean	53.26
Median	56
Min, Max	22, 73
Height (cm)	
п	50
Mean	164.3
Median	162
Min, Max	150.0, 182.0
Weight (kg)	
n	50
Mean	63.76
Median	63.5
Min, Max	42, 93
BMI (kg/m <sup>2</sup> )	
n	50
Mean	23.5
Median	23.34
Min, Max	16.41, 31.25
Medical history	
Esophagitis	48
Barrett's esophagus	1
Esophageal hiatal hernia	2
Superficial gastritis	8
Atrophic gastritis	2

EGD: Esophagogastroduodenoscopy; ds-MCE: Detachable string magnetically controlled capsule endoscopy; BMI: Body mass index.

86% (95%CI: 72.64%, 93.72%; Tables 3 and 4). Moreover, ds-MCE was also capable of clearly capturing and preserving images of the upper, middle, and lower esophagus and the dentate line with 100% integrity (Figure 3).

#### Assessment of the comfort and convenience of ds-MCE

The perceived comfort of the subjects who underwent EGD and ds-MCE examinations was assessed *via* a questionnaire. The overall comfort level of EGD was lower than that of the ds-MCE examination. Among them, 40 patients (80%) perceived that it was very comfortable or comfortable to undergo a ds-MCE examination, but only 19 patients (38.7%) perceived that it was very comfortable or comfortable to undergo an EGD. In addition, 7 patients (14.3%) considered undergoing EGD to be uncomfortable or very uncomfortable, and the data were significantly different between the two sets (P < 0.0001). Regarding the convenience of the ds-MCE examination, 25 patients (50%) found it very convenient, and no patients found it inconvenient (Figure 4 and Table 5).

Baishideng® WJG | https://www.wjgnet.com

Table 3 Detection endoscopy	n of esophageal lesions by esophagogastroduodenoscopy and detachable string magnetically controlled capsule				
		EGD			
		Positive	Negative	Total	
ds-MCE	Positive	18	4	22	
	Negative	3	25	28	
	Total	21	29	50	

EGD: Esophagogastroduodenoscopy; ds-MCE: Detachable string magnetically controlled capsule endoscopy.

Table 4 Detection of esophageal lesions by esophagogastroduodenoscopy and detachable string magnetically controlled capsule endoscopy

Primary endpoint	95%CI
Sensitivity	85.71% (62.64%, 96.24%)
Specificity	86.21% (67.43%, 95.49%)
Positive predictive value	81.82% (58.99%,94.01%)
Negative predictive value	89.29% (70.63%, 97.19%)
Diagnostic accuracy	86% (72.64%,93.72%)




Figure 3 Representative diseases observed on detachable string magnetically controlled capsule endoscopy and esophagogastroduodenoscopy. A-D: Detachable string magnetically controlled capsule endoscopy; E-H: Esophagogastroduodenoscopy.

## DISCUSSION

The evolution of diagnostic modalities in gastroenterology is driven by the dual requirements of improving patient comfort and expanding clinicians' diagnostic capabilities. In this context, the development of ds-MCE represents a detection innovation offering patients a less invasive option and improving compliance with recommended diagnostic procedures. Our clinical study evaluated the safety and effectiveness of ds-MCE, a novel capsule endoscopy technique, for diagnosing esophageal diseases. Our findings confirmed the clinical utility of ds-MCE as a feasible, accurate, and safe noninvasive diagnostic modality that may provide tremendous benefits in terms of patient comfort and diagnostic value.

Our study identified specific patient populations that may benefit from ds-MCE detection. These included individuals at lower risk for esophageal pathology who did not require tissue sampling, such as those undergoing routine surveillance for conditions such as Barrett's esophagus. Furthermore, ds-MCE presented a safer diagnostic option for

## Table 5 Perception analysis of subjects undergoing esophagogastroduodenoscopy and detachable string magnetically controlled capsule endoscopy examination

Projects EGD (n = 49) ds-MCE (n = 50)   Comfort level 0 Very comfort 6 13   1 Comfort 13 27	
0 Very comfort 6 13	
1 Comfort 13 27	
2 Tolerable 23 10	
3 Uncomfortable 3 0	
4 Very uncomfortable 4 0	
<i>P</i> value (Kruskal-Wallis test) < 0.0001	
Pain during the examination	
0 No 23 40	
1 A little 18 9	
2 Mild 7 1	
3 Serious 1 0	
4 Unbearable 0 0	
<i>P</i> value (Kruskal-Wallis test) 0.0004	
Discomfort during the examination	
0 No 14 31	
1 A little 21 14	
2 Mild 6 5	
3 Serious 7 0	
4 Unbearable 0 0	
<i>P</i> value (Kruskal-Wallis test) 0.0003	
Pain after the examination	
0 No 29 44	
1 A little 15 6	
2 Mild 5 0	
3 Serious 0 0	
4 Unbearable 0 0	
<i>P</i> value (Kruskal-Wallis test) 0.0008	
Discomfort after the examination	
0 No 23 42	
1 A little 18 7	
2 Mild 7 1	
3 Serious 1 0	
4 Unbearable 0 0	
<i>P</i> value (Kruskal-Wallis test) < 0.0001	
Convenience of inspection	
0 Very convenient 14 25	
1 Convenient 29 25	
2 Inconvenient 5 0	
3 Very inconvenient 1 0	



P value (Kruskal-Wallis test)

0.0074

EGD: Esophagogastroduodenoscopy; ds-MCE: Detachable string magnetically controlled capsule endoscopy.




Figure 4 Perception analysis of subjects undergoing esophagogastroduodenoscopy and detachable string magnetically controlled capsule endoscopy examination. A: Comfort level; B: Pain during the examination; C: Discomfort during the examination; D: Pain after the examination; E: Discomfort after the examination; F: Convenience of inspection. <sup>a</sup>P < 0.05, <sup>b</sup>P < 0.01, <sup>c</sup>P < 0.001, <sup>d</sup>P < 0.0001. ds-MCE: Detachable string magnetically controlled capsule endoscopy; EGD: Esophagogastroduodenoscopy.

patients at increased risk of complications from traditional endoscopy or sedation, including elderly individuals, pregnant women, and those with hemodynamic instability. This approach might meet the growing demand for patientfriendly diagnostics that minimize discomfort and risk.

The unique design of ds-MCE features a steerable capsule endoscope tethered by a detachable cord, which facilitates detailed examination of targeted esophageal areas. ds-MCE has demonstrated enhanced visualization rates across various esophageal segments, including the Z-line, with one study reporting significantly greater visualization than standard MCE[17]. This ability to obtain high-quality mucosal images without anesthesia represents a significant advance in patient-centered care, especially for those who are averse to traditional endoscopy, confirming the results of previous studies[14]. Comparative analysis with existing technologies such as the N-scope demonstrated that ds-MCE can enhance the patient experience by minimizing discomfort. However, direct comparative studies are needed to confirm these claims.

The design of ds-MCE allows for extended examination times within the esophagus, which is a substantial improvement over traditional capsule endoscopy and often provides limited esophageal imaging due to rapid transit. The limitations inherent in current methods underscore the need for innovation in diagnostic endoscopy. The ds-MCE system in our study addressed these limitations by allowing extended visualization of the esophageal mucosa and was particularly advantageous for assessing areas that are traditionally challenging to evaluate using standard endoscopic techniques

Despite these promising results, we acknowledge that our study has several limitations. First, patients with esophageal or gastric varices, which are commonly found in the lower esophagus, were not included, leaving it unclear whether their presence would alter the interpretation of the dentate line by ds-MCE. Second, most of the participants in this study were older adults with high levels of anxiety, and whether the comfort ratings would be similar for younger individuals needs to be investigated. Third, this pilot study recruited adults but did not address the potential challenges of capsule ingestion in groups known for high rates of swallowing difficulties, such as children. Fourth, ds-MCE takes longer than does EGD, with median examination times of 14.3 and 6.2 min, respectively [18]. Fifth, as a noninvasive procedure, ds-MCE does not provide histological samples or therapeutic options[19], and such considerations should be weighed within the broader clinical workflow and diagnostic strategy.

The clinical significance of false negatives in diagnosing esophageal pathology, especially for high-risk lesions such as early-stage malignancies, is profound. Although the sensitivity shown in our study is encouraging, the observed falsenegative rate of approximately 15% necessitates a cautious approach when considering ds-MCE for the detection of such critical conditions. Rigorous evaluation of the diagnostic accuracy of ds-MCE is required in subsequent research efforts.

Zaishideng® WJG | https://www.wjgnet.com

Prospective, multicenter studies with larger sample sizes and the inclusion of high-risk patient groups are essential for validating our findings. Such studies should also perform head-to-head comparisons with established gold-standard diagnostic procedures to firmly establish the reliability and clinical applicability of ds-MCE.

## CONCLUSION

In summary, our study highlights ds-MCE as a viable alternative to conventional EGD, as it has comparable efficacy in detecting esophageal lesions and has the potential to improve patient compliance in clinical settings. The continuous improvement and development of capsule endoscopy technology is expected to revolutionize the detection and management of upper gastrointestinal pathology [20]. To fully realize the potential of ds-MCE, future studies should include diverse patient groups, such as children and elderly individuals, to ensure robust and generalizable results. Investigating the cost-effectiveness of ds-MCE and its integration into standard diagnostic algorithms will be pivotal for its widespread adoption in clinical practice. Preliminary cost analysis models, accounting for potential reductions in anesthesia use and increased patient throughput, can offer valuable insights into the economic viability of ds-MCE. Through these coordinated research efforts, ds-MCE is likely to redefine the paradigm of esophageal disease diagnosis in the coming years.

## ARTICLE HIGHLIGHTS

#### Research background

Esophagogastroduodenoscopy (EGD), the gold standard for diagnosing esophageal diseases, is invasive and can cause patient discomfort. Magnetically controlled capsule endoscopy (MCE) is a noninvasive alternative, yet its application for esophageal conditions is limited.

#### Research motivation

This study was motivated by the need to enhance patient compliance and comfort during esophageal examinations while maintaining high diagnostic accuracy.

#### Research objectives

The primary aim was to assess the safety and diagnostic efficacy of detachable string MCE (ds-MCE) for identifying esophageal diseases.

#### Research methods

A prospective, blinded, self-controlled clinical study was conducted to compare ds-MCE with EGD across various diagnostic parameters and patient tolerance levels in 50 subjects with known esophageal diseases.

#### Research results

Ds-MCE demonstrated diagnostic accuracy comparable to that of EGD (86%) and was preferred by patients, with 80% reporting a comfortable or very comfortable experience.

#### Research conclusions

Ds-MCE is an effective and patient-friendly diagnostic tool for esophageal diseases, offering a noninvasive alternative to traditional EGD with similar diagnostic outcomes.

#### Research perspectives

Further studies should expand patient demographics to validate the utility of ds-MCE across various populations and evaluate its cost-effectiveness for broader clinical adoption.

## FOOTNOTES

Co-corresponding authors: Jin-Shui Zhu and Jing Zhang.

Author contributions: Zhang J, Zhu JS and Da W conceptualized and designed the research; Da W, Yang YL and Qin HW screened patients and acquired clinical data; Chen ZY and Fan HN collected the clinical data; Yu Y performed Data analysis; Yang YL wrote the paper. All the authors have read and approved the final manuscript. Both Zhang J and Zhu JS have played important and indispensable roles in the experimental design, data interpretation and manuscript preparation as the co-corresponding authors. Zhu JS applied for and obtained the funds for this research project. Zhu JS conceptualized, designed, and supervised the whole process of the project. Zhu JS searched the literature, revised the early version of the manuscript with the focus on the safety of detachable string magnetically controlled capsule endoscopy (ds-MCE) for identifying esophageal diseases. Zhang J was instrumental and responsible for data reanalysis and re-interpretation, figure plotting, comprehensive literature search, preparation and submission of the current version of the manuscript with a new focus on the efficacy of ds-MCE for the diagnosis of esophageal diseases. This collaboration between Zhang J and



Yang YL et al. A prospective, blind clinical study

Zhu JS is crucial for the publication of this manuscript and other manuscripts still in preparation.

Supported by the Science and Technology Commission of Shanghai, No. 18DZ1930309.

Institutional review board statement: The study was reviewed and approved by the Shanghai Sixth People's Hospital Institutional Review Board, No. 2019-082-(1).

Clinical trial registration statement: This study was registered with Clinical Trials.gov. The registration identification number is NCT05469152.

Informed consent statement: All subjects understood and agreed to the study protocol and voluntarily signed the informed consent form

**Conflict-of-interest statement:** All the authors report no relevant conflicts of interest for this article.

Data sharing statement: No additional data are available.

CONSORT 2010 statement: The authors have read the CONSORT 2010 statement, and the manuscript was prepared and revised according to the CONSORT 2010 statement.

**Open-Access:** This article is an open-access article that was selected by an in-house editor and fully peer-reviewed by external reviewers. It is distributed in accordance with the Creative Commons Attribution NonCommercial (CC BY-NC 4.0) license, which permits others to distribute, remix, adapt, build upon this work non-commercially, and license their derivative works on different terms, provided the original work is properly cited and the use is non-commercial. See: https://creativecommons.org/Licenses/by-nc/4.0/

#### Country/Territory of origin: China

**ORCID number:** Jin-Shui Zhu 0000-0002-3175-568X; Jing Zhang 0000-0002-9412-3567.

S-Editor: Li L L-Editor: A P-Editor: Chen YX

## REFERENCES

- Zagari RM, Fuccio L, Wallander MA, Johansson S, Fiocca R, Casanova S, Farahmand BY, Winchester CC, Roda E, Bazzoli F. Gastrooesophageal reflux symptoms, oesophagitis and Barrett's oesophagus in the general population: the Loiano-Monghidoro study. Gut 2008; 57: 1354-1359 [PMID: 18424568 DOI: 10.1136/gut.2007.145177]
- Eusebi LH, Cirota GG, Zagari RM, Ford AC. Global prevalence of Barrett's oesophagus and oesophageal cancer in individuals with gastro-2 oesophageal reflux: a systematic review and meta-analysis. Gut 2021; 70: 456-463 [PMID: 32732370 DOI: 10.1136/gutjnl-2020-321365]
- Shaheen O, Ghibour A, Alsaid B. Esophageal Cancer Metastases to Unexpected Sites: A Systematic Review. Gastroenterol Res Pract 2017; 3 2017: 1657310 [PMID: 28659974 DOI: 10.1155/2017/1657310]
- 4 Thrift AP. Global burden and epidemiology of Barrett oesophagus and oesophageal cancer. Nat Rev Gastroenterol Hepatol 2021; 18: 432-443 [PMID: 33603224 DOI: 10.1038/s41575-021-00419-3]
- Cardey J, Le Gall C, Michaud L, Dabadie A, Talbotec C, Bellaiche M, Lamireau T, Mas E, Bridoux-Henno L, Viala J, Restier-Miron L, 5 Lachaux A. Screening of esophageal varices in children using esophageal capsule endoscopy: a multicenter prospective study. Endoscopy 2019; 51: 10-17 [PMID: 30184608 DOI: 10.1055/a-0647-1709]
- Bouchard S, Ibrahim M, Van Gossum A. Video capsule endoscopy: perspectives of a revolutionary technique. World J Gastroenterol 2014; 6 20: 17330-17344 [PMID: 25516644 DOI: 10.3748/wjg.v20.i46.17330]
- Eliakim R, Sharma VK, Yassin K, Adler SN, Jacob H, Cave DR, Sachdev R, Mitty RD, Hartmann D, Schilling D, Riemann JF, Bar-Meir S, 7 Bardan E, Fennerty B, Eisen G, Faigel D, Lewis BS, Fleischer DE. A prospective study of the diagnostic accuracy of PillCam ESO esophageal capsule endoscopy vs conventional upper endoscopy in patients with chronic gastroesophageal reflux diseases. J Clin Gastroenterol 2005; 39: 572-578 [PMID: 16000923 DOI: 10.1097/01.mcg.0000170764.29202.24]
- Gralnek IM, Adler SN, Yassin K, Koslowsky B, Metzger Y, Eliakim R. Detecting esophageal disease with second-generation capsule 8 endoscopy: initial evaluation of the PillCam ESO 2. Endoscopy 2008; 40: 275-279 [PMID: 18389444 DOI: 10.1055/s-2007-995645]
- 9 Hosoe N, Naganuma M, Ogata H. Current status of capsule endoscopy through a whole digestive tract. Dig Endosc 2015; 27: 205-215 [PMID: 25208463 DOI: 10.1111/den.12380]
- Duvvuri A, Desai M, Vennelaganti S, Higbee A, Gorrepati VS, Dasari C, Chandrasekar VT, Vennalaganti P, Kohli D, Sathyamurthy A, Rai T, 10 Sharma P. Diagnostic accuracy of a novel third generation esophageal capsule as a non-invasive detection method for Barrett's esophagus: A pilot study. J Gastroenterol Hepatol 2021; 36: 1222-1225 [PMID: 32996655 DOI: 10.1111/jgh.15283]
- Park J, Cho YK, Kim JH. Current and Future Use of Esophageal Capsule Endoscopy. Clin Endosc 2018; 51: 317-322 [PMID: 30078304 DOI: 11 10.5946/ce.2018.101]
- Zou WB, Hou XH, Xin L, Liu J, Bo LM, Yu GY, Liao Z, Li ZS. Magnetic-controlled capsule endoscopy vs. gastroscopy for gastric diseases: a 12 two-center self-controlled comparative trial. Endoscopy 2015; 47: 525-528 [PMID: 25590177 DOI: 10.1055/s-0034-1391123]
- 13 Liao Z, Hou X, Lin-Hu EQ, Sheng JQ, Ge ZZ, Jiang B, Hou XH, Liu JY, Li Z, Huang QY, Zhao XJ, Li N, Gao YJ, Zhang Y, Zhou JQ, Wang XY, Liu J, Xie XP, Yang CM, Liu HL, Sun XT, Zou WB, Li ZS. Accuracy of Magnetically Controlled Capsule Endoscopy, Compared With



WJG | https://www.wjgnet.com

Conventional Gastroscopy, in Detection of Gastric Diseases. Clin Gastroenterol Hepatol 2016; 14: 1266-1273.e1 [PMID: 27211503 DOI: 10.1016/j.cgh.2016.05.013]

- Beg S, Card T, Warburton S, Rahman I, Wilkes E, White J, Ragunath K. Diagnosis of Barrett's esophagus and esophageal varices using a 14 magnetically assisted capsule endoscopy system. Gastrointest Endosc 2020; 91: 773-781.e1 [PMID: 31678203 DOI: 10.1016/j.gie.2019.10.031]
- Chen YZ, Pan J, Luo YY, Jiang X, Zou WB, Qian YY, Zhou W, Liu X, Li ZS, Liao Z. Detachable string magnetically controlled capsule 15 endoscopy for complete viewing of the esophagus and stomach. Endoscopy 2019; 51: 360-364 [PMID: 30822803 DOI: 10.1055/a-0856-6845]
- Chen XY, Da W, Liang R, Fan HN, Yi YC, Chen M, Qin HW, Zhang J, Zhu JS. The Detective Value of Magnetically Controlled Robotic 16 Capsule Endoscopy in Patients With Suspected Small Intestinal Disease. Front Med (Lausanne) 2021; 8: 610563 [PMID: 34113626 DOI: 10.3389/fmed.2021.610563]
- Song J, Bai T, Zhang L, Xiang XL, Xie XP, Hou XH. Better view by detachable string magnetically controlled capsule endoscopy for 17 esophageal observation: a retrospective comparative study. Dis Esophagus 2020; 33 [PMID: 32008045 DOI: 10.1093/dote/doz104]
- Jiang X, Qian YY, Liu X, Pan J, Zou WB, Zhou W, Luo YY, Chen YZ, Li ZS, Liao Z. Impact of magnetic steering on gastric transit time of a 18 capsule endoscopy (with video). Gastrointest Endosc 2018; 88: 746-754 [PMID: 30005825 DOI: 10.1016/j.gie.2018.06.031]
- Wang S, Huang Y, Hu W, Mao H, McAlindon ME, Liu Y, Yang L, Zhang C, Xu M, He C, Dang T, Wu B, Ji D, Zhang L, Mao X, Liu C, Xu 19 D, Li Y, Li G, Han J, Lv F, Liang X, Jin S, Zhang S, Tai FWD, Xu Q, Yang C, Wang G, Wang L, Li B, Yang H, Xie P, Deng L, Ren L, Chang Z, Wang X, Wang S, Gao X, Li J, Zhu L, Wang F, Zhang G, Jiang X, Pan J, Meng W, Li X, Hou J, Dray X, Liao Z, Qi X. Detachable string magnetically controlled capsule endoscopy for detecting high-risk varices in compensated advanced chronic liver disease (CHESS1801): A prospective multicenter study. Lancet Reg Health West Pac 2021; 6: 100072 [PMID: 34327406 DOI: 10.1016/j.lanwpc.2020.100072]
- Alsayid M, Melson J. Will magnet-assisted capsule endoscopy become a viable screening tool for Barrett's esophagus and esophageal varices? 20 Gastrointest Endosc 2020; 91: 782-784 [PMID: 32204813 DOI: 10.1016/j.gie.2019.12.015]





## Published by Baishideng Publishing Group Inc 7041 Koll Center Parkway, Suite 160, Pleasanton, CA 94566, USA Telephone: +1-925-3991568 E-mail: office@baishideng.com Help Desk: https://www.f6publishing.com/helpdesk https://www.wjgnet.com

